



ENJOY CLAIRE

cosmetic device for the treatment of cellulite

model 3.0



user manual

Contents

Introduction.....	2
Properties of the device.....	3
Device components	4
Description of control panel	5
Description of rear panel	6
Instruction on use	7
Maintenance	10
IMPORTANT INSTRUCTION	10
Technical data, failures	12

Introduction

ENJOY Claire is a magnetotherapeutic device that offers ladies a comfortable way of getting rid of the unpleasant “orange peel” appearance. Its therapeutically acting components are **planar pulsed magnetic field and heat**, which accelerate the metabolism of fatty cells and remove waste products and toxins via the lymphatic system. The two types of energy are generated by a planar applicator and penetrate into the body only to a depth of a few millimeters.

It is sufficient to remain in sitting or lying position and let the pleasant heat of the applicator act for at least **an hour a day regularly for a period of 20 days**. To facilitate the lymph draining through the lymphatic system it is recommended to have the legs slightly raised during application. **A pronounced change can only be expected after a fortnight**. For a prolonged preservation of the result achieved, it is recommended to repeat the treatment in shorter or longer intervals, as may be necessary. The same as in skin cosmetics, in the case of cellulite too, it is necessary to avoid blemishes by preventive stimulation of metabolism, which this device does.

Properties of the device

- adjustable magnetic field intensity (MFI) on the applicator surface
- control of 43 °C max. temp. as a protection against skin overheating
- automatic switching off after 60 minutes' operation or in case of failure
- using a very low voltage in order to increase user safety
- high effectiveness, which reduces operating costs and environmental damage
- possibility of washing the applicator cover in current detergents


Device components

1. ENJOY Claire control unit, inclusive of the applicator
2. Power supply 24 V/2.7 A SYS 1183-6524
3. Supply cable
4. Textile cover of applicator
5. Disposable absorbent pad
6. User manual
7. Carrier bag




Description of control panel

- a. MFI degree 1 on/off (applicator surface temperature ca. 39 °C)
- b. MFI degree 2 on/off (applicator surface temperature ca. 41 °C)
- c. MFI degree 3 on/off (applicator surface temperature ca. 43 °C)
- d. Signal lamp of IMF degree 1
- e. Signal lamp of IMF degree 2
- f. Signal lamp of IMF degree 3



Description of rear panel

- A. Main switch
- B. Power supply connector
- C. Program facility connector
- D. Applicator output



Instruction on use

CONNECTING AND PUTTING INTO OPERATION

1. Place the control unit (1) on a firm base.
2. Put the textile cover (4) on the applicator and close the zippers, the inscription "APPLICATION SIDE" being on the outside.
3. The textile cover does not protect the applicator against dampness (sweat); apply therefore a disposable absorbent pad (5), which is placed between the textile cover and the body.
4. Place the applicator on the therapy site (a chair, armchair, bed, etc.). Take care that during application the inscription "APPLICATION SIDE" is towards the body.
5. Connect the power supply (2) to the connector (B) of the control unit (1). Do not cover the power supply lest it should be damaged.
6. Connect the supply cable (3) to the power supply (2) and then to a mains socket (technical data).

To speed up the healing effect, we recommend switching on and preheating the device some 30 minutes prior to application. Cover the applicator with a blanket to prevent its cooling. After the above period switch the device off and proceed as given under 2 and 3 below.

PREPARATION FOR THERAPY

1. Sit down on the applicator and tighten the belts with Velcro fasteners around the thighs. Be aware that any layer of clothes reduces the healing effect.
2. During application, lift your legs to a slightly higher position.

SWITCHING THE DEVICE ON AND THERAPY PROPER

1. Turn the main switch (A) to position **I** on the rear panel of the control unit (1) and check whether all the free control/signal lamps/lights (d, e, f) on the control panel are blinking.
2. Choose one of the MFI degrees by pushing the respective button (a, b, c) on the control panel (1). From now on, the device is in operation and the signal lamp of the chosen MFI degree must give continuous light.
3. After an hour, the therapy will automatically be terminated; you will be warned by a sound signal. You can terminate the therapy before due time by repeatedly pushing the same button.
4. During the therapy you can change the MFI degree to your liking.

SWITCHING THE DEVICE OFF

1. When the therapy is over, turn the switch (A) to position 0 and check whether all the three signal lamps (d, e and f) go out.
2. Withdraw the supply cable (3) from the mains socket.

Maintenance

- Before cleaning or washing the textile cover take it off the applicator. **Set the washing temperature to 30 °C at the most.** After cleaning or washing **dry the cover thoroughly but do not iron it.**
- Do not wash the applicator or the control unit.
- Do not use aggressive chemicals when cleaning or washing the cover.

IMPORTANT INSTRUCTIONS

1. Read the user manual before using the device for the first time. Keep the user manual for later reference.
2. The device does not substitute any therapeutic facility designed for medical application in hospitals.
3. If you are in doubt as to your health condition, contact your doctor before using the device. Do not use the device during pregnancy.
4. Do not use the device in your sleep and do not apply it to people while they are unconscious.
5. Do not use the device for longer than 2 hours without a break of several hours.
6. Do not apply the applicator to skins that are oversensitive, diseased or under treatment. Terminate the application immediately you notice any unpleasant sensation. During application, check occasionally the condition of the skin exposed to treatment.
7. Do not open or modify any part of the device.
8. Do not disconnect the device by pulling the mains cable.

9. Use only original accessories supplied by the manufacturers, the power supply in particular.
10. Do not expose the device to direct sunshine, high temperatures or chemicals.

CAUTION:

The device is not designed for use by persons (children included) whose physical and mental disabilities or lack of experience and knowledge prevents them from using the device safely unless they are supervised or have been advised about the functioning of the device by a person responsible for their safety. The device should be kept out of the reach of children. The device is not designed for playing with.

Technical data, failures

ENJOY CLAIRE

Model:	Char-M05N_V3
Power supply:	DC 24 V/2.7 A
Power input:	65 W
Therapy duration setting:	60 min
Maximum applicator temperature:	43 °C
Trans dimensions (l/w/h):	(450/450/110) mm

POWER SUPPLY UNIT

Input voltage:	AC (110-240) V, (50-60) Hz
Output voltage:	DC 24 V
Output current:	2.7 A
Output power:	65 W
Protection class:	II, double insulation
Output connector:	(2,1 x 5,5) mm

Failures

If after switching on the device the three signal lamps (d, e and f) are not blinking or if the chosen required MFI degree cannot be obtained and a continuous sound signal is heard or if the device gets disconnected prematurely, switch the device off by the main switch (A) and disconnect it from the mains. After five seconds, connect the device and switch it on again. If the problem persists, contact the vendor.



ENJOY spol. s r.o.

Loosova 355/12, 638 00 Brno, ID: 46968270, Czech Republic

www.enjoy.cz info@enjoy.cz